

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	3D
Patient Position	Supine	Pulse Sequence	SPGR
Coil Configuration	8HRBRAIN	Imaging Options	EDR, Fast, IrP, Asset
Plane	SAGITTAL		
Series Description	FSPGR_Anat		
SCAN TIMING		SCAN RANGE	
Flip Angle	8	FOV	25.6
TE	Min Full	Slice Thickness	1.0
Number of Echoes	1	Location per Slab	176
TI	600	Overlap Locations	0
Receiver Bandwidth	31.25		
IMAGE ENHANCE		ACQ TIMING	
Filter Choice	None	Freq	256
		Phase	192
		Freq DIR	S/I
		NEX	1.00
		Phase FOV	1.00
		Auto Shim	Auto
		Phase Correction	No
GATING/TRIGGER		USER CVS	
Auto Trigger Type	Off	User CV23	100.00
FMRI		MULTI-PHASE	
Initial State	Control	# of Acquisition	0
PSD Trigger	Internal	Seperate Series	0
Slice Order	Interleaved	Mask Phase	0
View Order	Bottom/Up	Mask Pause	0
# of Repetitions REST	0		
# of Repetitions ACTIVE	0		
SAT		TRICKS	
Tag Type	None	Pause On/Off	On
		Auto Subtract	0
		Auto SCIC	Off
CONTRAST			
Contrast Yes/No	No		

FSPGR_Anat

FSPGR_Anat