

Resting_State	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Imaging Mode	2D
	Patient Position	Supine	Pulse Sequence	Gradient Echo
	Coil Configuration	32 Ch Head Coil	Imaging Options	EPI, Asset, FMRI
	Plane	AXIAL	n/a	9
	Series Description	Resting_State	SCAN RANGE	
	SCAN TIMING		FOV	22.0
	Flip Angle	75	Slice Thickness	3.4
	TE	30.0	Slice Spacing	0.0
	TR	2000.0	End Location 2	0.0
	Number of Shots	1	End Location 3	0.0
	IMAGE ENHANCE		Center of Location Start	0.0
	Filter Choice	None	Center of Location End	0.0
	GATING/TRIGGER		Number of Slices	42
	Heart Beat per Minute Mode	0	Slice for 3 Plane Localizaion	0
	Auto Trigger Type	Off	Slice for 3 Plane Localizaion	0
	Auto Trigger Window	Off	Slice for 3 Plane Localizaion	0
	FMRI		Slice for 3 Plane Localizaion	0
	Brain Wave Real Time	2	Space per Plane 1	0.0
	Paradigm String	FACE_LOC	Space per Plane 2	0.0
	Paradigm UID	1.2.840.113819.3.69136390 610.1372432428.356	Table Delta	0.00
	Initial State	Control	ACQ TIMING	
	PSD Trigger	Internal	Freq	64
	Slice Order	Interleaved	Phase	64
	View Order	Bottom/Up	Freq DIR	R/L
	# of Repetitions REST	6	NEX	1.00
	# of Repetitions ACTIVE	6	Phase FOV	1.00
	# of Dummy Acquisition	5	Auto Shim	Auto
SAT		Phase Correction	Yes	
Tag Type	None	USER CVS		
TRICKS		User CV0	1.00	
Pause On/Off	On	User CV Mask	1	
Auto Subtract	0	MULTI-PHASE		
Auto SCIC	Off	Slice per Location	295	
		Delay after Acquisition	0	
		Seperate Series	0	
		Delay after Acquisition without AV	0	
		Mask Phase	0	
		Mask Pause	0	
		DIFFUSION		
		Recon All Images	On	
		CONTRAST		
		Contrast Yes/No	No	
		Contrast Amount	Yes	

Resting_State