

INFO PAGE		GEOMETRY		CONTRAST	
Total scan duration	08:08.3	Nucleus	H1	Scan type	Imaging
Rel. signal level (%)	100	Coil selection	SENSE-Head-8	Scan mode	3D
Act. TR/TE (ms)	8.0 / 3.7	element selection	SENSE	technique	FFE
ACQ matrix M x P	256 x 200	connection	d	Contrast enhancement	T1
ACQ voxel MPS (mm)	1.00 / 1.00 / 1.00	Dual coil	no	Acquisition mode	cartesian
REC voxel MPS (mm)	1.00 / 1.00 / 1.00	CLEAR	yes	Fast Imaging mode	TFE
Scan percentage (%)	100	body tuned	yes	3D non-selective	no
TFE shots	140	FOV FH (mm)	256	shot mode	multishot
TFE dur. shot / acq (ms)	1814.0 / 1610.0	RL (mm)	200	TFE factor	200
Min. TI delay	843.2495	AP (mm)	200	3D free factor	no
Act. WFS (pix) / BW (Hz)	2.269 / 191.5	Voxel size FH (mm)	1	startup echoes	default
Min. WFS (pix) / Max. BW (Hz)	0.478 / 909.7	RL (mm)	1	shot interval	user defined
SAR / head	< 5 % / 0.2 W/kg	AP (mm)	1	(ms)	3500
Whole body / level	0.0 W/kg / normal	Recon voxel size (mm)	1	profile order	linear
B1 rms [uT]	0.5456858	RFOV (%)	78.125	turbo direction	Y
PNS / level	60 % / normal	Fold-over suppression	no	Echoes	1
Sound Pressure Level (dB)	12.25113	Slice oversampling	default	partial echo	no
MOTION		Matrix scan	256	shifted echo	no
Cardiac synchronization	no	reconstruction	256	TE	shortest
Respiratory compensation	no	Scan percentage (%)	100	Flip angle (deg)	8
Navigator respiratory comp	no	SENSE	yes	TR	shortest
Flow compensation	no	P reduction (RL)	1	Halfscan	no
fMRI echo stabilisation	no	P os factor	1	Water-fat shift	user defined
Motion smoothing	no	S reduction (AP)	2	(pixels)	3.5
NSA	1	k-t BLAST	no	Shim	auto
DYN/ANG		Overcontiguous slices	no	Fat suppression	no
Angio / Contrast enh.	no	Stacks	1	Water suppression	no
Quantitative flow	no	slices	200	TFE prepulse	invert
CENTRA	no	slice thickness (mm)	1	slice selection	no
Manual start	no	slice orientation	coronal	delay	user defined
Dynamic study	no	fold-over direction	RL	To Exc (ms)	1000
Arterial Spin labeling	no	fat shift direction	F	PSIR	no
vascular crushing	no	Chunks	1	tfe switch	no
Mag reset	no	PlanAlign	no	MAGIC enable	no
POST/PROC		REST slabs	0	MTC	no
Preparation phases	auto	Interactive positioning	no	T2prep	no
Power optimization	auto	OFFC/ANG		Research prepulse	no
Manual Offset Freq.	no	Stacks	1	Diffusion mode	no
SENSE ref. scan	no	Stack Offc. AP (P=+mm)	-20.90301	SAR mode	high
B0 field map/Dixon	no	RL (L=+mm)	0	B1 mode	default
B1 field map	no	FH (H=+mm)	15.46823	PNS mode	low
MIP/MPR	no	Ang. AP (deg)	0	Gradient mode	maximum
Images	M, no, no, no	RL (deg)	6.465102	SofTone mode	no
Autoview image	M	FH (deg)	0		
Calculated images	no, no, no, no				
Reference tissue	Grey matter				
Preset window contrast	soft				
Reconstruction mode	immediate				
Save raw data	no				
Hardcopy protocol	no				
Ringing filtering	rectangular				
Geometry correction	default				
Elliptical k-space shutter	default				