

# SIEMENS MAGNETOM Allegra syngo MR A30

\\USER\Dr. O'Hearn\Faces\Faces\TC\_t1\_mprage\_sag\_ns

Scan Time: 8:59    Voxel size: 1.1x1.1x1.1 [mm]    Rel. SNR: 1.00    SIEMENS: tfl

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 [%]
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0 [deg]
Phase oversampling	0 [%]
Slice oversampling	0 [%]
Slices per slab	176
FoV read	269 [mm]
FoV phase	100.0 [%]
Slice thickness	1.05 [mm]
TR	2100 [ms]
TE	3.93 [ms]
Averages	1
Concatenations	1
Filter	Elliptical filter
Coil elements	HE

## Contrast

Magn. preparation	Non-sel. IR
TI	1000 [ms]
Flip angle	7 [deg]
Fat suppr.	None
Water suppr.	None
-----	
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1

## Resolution

Base resolution	256
Phase resolution	100 [%]
Slice resolution	100 [%]
Phase partial Fourier	Off
Slice partial Fourier	Off
Filter 1	
Raw filter	Off
Filter 2	
Large FoV	Off
Filter 3	
Normalize	Off
Filter 4	
Elliptical filter	On
Filter 5	
Image Filter	Off
Trajectory	Cartesian
Interpolation	0
-----	
PAT mode	None

## Geometry

Multi-slice mode	Single shot
Series	Ascending

## System

Save uncombined	0
Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H

## Head 3T / HE

Head 3T / HE	1
-----	
Shim mode	Standard
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	Isocenter
Orientation	Sagittal
Rotation	0 [deg]
F >> H	269 [mm]
A >> P	269 [mm]
R >> L	185 [mm]

## Physio

1st Signal/Mode	None
-----	
Dark blood	0
-----	
Resp. control	Off

## Inline

Subtract	0
Std-Dev-Sag	0
Std-Dev-Cor	0
Std-Dev-Tra	0
Std-Dev-Time	0
MIP-Sag	0
MIP-Cor	0
MIP-Tra	0
MIP-Time	0
Save original images	1

## Sequence

Introduction	1
Dimension	3D
Elliptical scanning	0
Asymmetric echo	Allowed
Bandwidth	130 [Hz/Px]
Echo spacing	9.4 [ms]
-----	
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	1