

SIEMENS MAGNETOM Allegra syngo MR 2004A

\USER\INVESTIGATORS\Castellanos\Child-AdultR01\dti_64dir_3x3x3

+ Scan Time: 5:43 Voxel size: 3.0x3.0x3.0 [mm] Rel. SNR: 1.00 USER: cbi_ep2d_diff_1_4

Routine	Physio	
Slice group 1	1st Signal/Mode	None
Slices	Diffusion mode	MDDW
Dist. factor	Diff. weightings	2
Position	b-value[1]	0 [s/mm ²]
Orientation	b-value[2]	1000 [s/mm ²]
Phase enc. dir.	Diff. weighted images	1
Rotation	Trace weighted images	0
Phase oversampling	Average ADC maps	0
FoV read	Individual ADC maps	0
FoV phase	Noise level	40
Slice thickness	Diff. directions	64
TR	Sequence	
TE	Introduction	0
Averages	Averaging mode	Long term
Concatenations	Bandwidth	3720 [Hz/Px]
Filter	Free echo spacing	0
Coil elements	Echo spacing	0.32 [ms]
Contrast	EPI factor	64
MTC	RF pulse type	Normal
Magn. preparation	Gradient mode	Fast
Reconstruction	Image Reconstruction	Save Raw
Fat suppr.	Field Map Mode	FMap Off
Measurements	Readout Direction	Normal
Delay in TR		
Resolution		
Base resolution		
Phase resolution		
Phase partial Fourier		
Filter 1		
Raw filter		
Filter 2		
Large FoV		
Filter 3		
Normalize		
Filter 4		
Elliptical filter		
Interpolation		
Geometry		
Multi-slice mode	Interleaved	
Series	Interleaved	
Special sat.	None	
System		
Scan at current TP	1	
MSMA	S - C - T	
Sagittal	R >> L	
Coronal	A >> P	
Transversal	F >> H	
Head 3T / HE	1	
Shim mode	Standard	
Confirm freq. adjustment	0	
Assume Silicone	0	
Ref. amplitude [1H]	140.000 [V]	
Adjust volume		
Position	Isocenter	
Orientation	Transversal	
Rotation	90 [deg]	
A >> P	192 [mm]	
R >> L	192 [mm]	
F >> H	150 [mm]	