

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	3D
Patient Position	Supine	Pulse Sequence	SPGR
Coil Configuration	32 Ch Head Coil	Imaging Options	EDR, Fast, IrP, Asset
Plane	SAGITTAL	n/a	6
Series Description	Sag_T1_FSPGR_BRAVO	<b>SCAN RANGE</b>	
SCAN TIMING		FOV	25.6
Flip Angle	12	Slice Thickness	1.0
TE	Min Full	Location per Slab	180
Number of Echoes	1	Overlap Locations	0
TI	450	GRXOPT	0
Receiver Bandwidth	31.25	Start Location 1	L85.8
IMAGE ENHANCE		End Location 1	R89.2
Filter Choice	None	End Location 2	A27.6
GATING/TRIGGER		End Location 3	S3.5
Heart Beat per Minute Mode	0	Center of Location Start	A27.6
Auto Trigger Type	Off	Center of Location End	S3.5
Auto Trigger Window	Off	Number of Slices	1
FMRI		Slice for 3 Plane Localizaion	0
PSD Trigger	Internal	Slice for 3 Plane Localizaion	0
Slice Order	Interleaved	Slice for 3 Plane Localizaion	0
View Order	Bottom/Up	Slice for 3 Plane Localizaion	0
# of Repetitions REST	0	Space per Plane 1	0.0
# of Repetitions ACTIVE	0	Space per Plane 2	0.0
SAT		Table Delta	0.00
Tag Type	None	ACQ TIMING	
TRICKS		Freq	256
Pause On/Off	On	Phase	256
Auto Subtract	0	Freq DIR	S/I
Auto SCIC	Off	NEX	1.00
		Phase FOV	1.00
		Auto Shim	Auto
		Phase Correction	No
		USER CVS	
		User CV23	100.00
		In-range AutoTR Control	0
		User CV Mask	8388688
		MULTI-PHASE	
		Seperate Series	0
		Mask Phase	0
		Mask Pause	0
		DIFFUSION	
		Recon All Images	On
		CONTRAST	
		Contrast Yes/No	No
		Contrast Amount	Yes

Sag\_T1\_FSPGR\_BRAVO

Sag\_T1\_FSPGR\_BRAVO