

SIEMENS MAGNETOM TrioTim syngo MR B15

\USER\BNU_cst\Zang.Yufeng\peng.zulai_6thHos_20090805\DTI_64_18x18x25_p2_opt

TA: 8:11 PAT: 2 Voxel size: 1.8x1.8x2.5 mm Rel. SNR: 1.00 SIEMENS: ep2d_diff

| Properties | | Special sat. | None |
|---|------------------|--------------------------|------------------------|
| Prio Recon | Off | | |
| System | | | |
| Before measurement | | Body | Off |
| After measurement | | HEP | On |
| Load to viewer | On | HEA | On |
| Inline movie | Off | | |
| Auto store images | On | | |
| Load to stamp segments | Off | Positioning mode | FIX |
| Load images to graphic segments | Off | Table position | H |
| Auto open inline display | Off | Table position | 0 mm |
| AutoAlign Spine | Off | MSMA | S - C - T |
| Start measurement without further preparation | On | Sagittal | R >> L |
| Wait for user to start | On | Coronal | A >> P |
| Start measurements | single | Transversal | F >> H |
| Routine | | Coil Combine Mode | Adaptive Combine |
| Slice group 1 | | Auto Coil Select | Default |
| Slices | 49 | | |
| Dist. factor | 0 % | Shim mode | Standard |
| Position | R1.4 A23.1 H10.8 | Adjust with body coil | Off |
| Orientation | Transversal | Confirm freq. adjustment | Off |
| Phase enc. dir. | A >> P | Assume Silicone | Off |
| Rotation | 0.00 deg | Ref. amplitude 1H | 335.349 V |
| Phase oversampling | 0 % | Adjustment Tolerance | Auto |
| FoV read | 230 mm | Adjust volume | |
| FoV phase | 100.0 % | Position | R1.4 A23.1 H10.8 |
| Slice thickness | 2.5 mm | Orientation | Transversal |
| TR | 7200 ms | Rotation | 0.00 deg |
| TE | 104 ms | R >> L | 230 mm |
| Averages | 1 | A >> P | 230 mm |
| Concatenations | 1 | F >> H | 123 mm |
| Filter | None | | |
| Coil elements | HEA;HEP | Physio | |
| | | 1st Signal/Mode | None |
| | | Resp. control | Off |
| Contrast | | Diff | |
| MTC | Off | Diffusion mode | MDDW |
| Magn. preparation | None | Diff. weightings | 2 |
| Fat suppr. | Fat sat. | b-value 1 | 0 s/mm ² |
| | | b-value 2 | 1000 s/mm ² |
| Averaging mode | Long term | Diff. weighted images | On |
| Reconstruction | Magnitude | Trace weighted images | Off |
| Delay in TR | 0 ms | Average ADC maps | On |
| Multiple series | Off | Individual ADC maps | Off |
| Resolution | | FA maps | On |
| Base resolution | 128 | Mosaic | On |
| Phase resolution | 100 % | Tensor | On |
| Phase partial Fourier | Off | Noise level | 40 |
| Interpolation | Off | Diff. directions | 64 |
| | | Sequence | |
| PAT mode | GRAPPA | Introduction | On |
| Accel. factor PE | 2 | Bandwidth | 1396 Hz/Px |
| Ref. lines PE | 39 | Free echo spacing | Off |
| Matrix Coil Mode | Auto (Triple) | Echo spacing | 0.82 ms |
| Reference scan mode | Separate | | |
| Distortion Corr. | Off | EPI factor | 128 |
| Prescan Normalize | Off | RF pulse type | Normal |
| Raw filter | On | Gradient mode | Fast* |
| Elliptical filter | Off | | |
| Hamming | Off | | |
| Geometry | | | |
| Multi-slice mode | Interleaved | | |
| Series | Interleaved | | |