

Protocol: EPILEPSY SETTINGS

3D AXIAL

PATIENT POSITION	
Patient Entry	Head First
Patient Position	Supine
Coil Configuration	8HRBRAIN
Plane	AXIAL
Series Description	3D AXIAL
SCAN TIMING	
Flip Angle	8
TE	Min Full
Number of Echoes	1
TI	650
Receiver Bandwidth	31.25
IMAGE ENHANCE	
Filter Choice	None
GATING/TRIGGER	
Auto Trigger Type	Off
FMRI	
PSD Trigger	Internal
Slice Order	Interleaved
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0
SAT	
Tag Type	None
TRICKS	
Pause On/Off	On
Auto Subtract	0
Auto SCIC	Off
MULTI-STATION	
Number of Stations	0
Station Number	0

3D AXIAL

IMAGING PARAMETERS	
Gradient Mode	ZOOM
Imaging Mode	3D
Pulse Sequence	SPGR
Imaging Options	EDR, Fast, IrP
Acceleration Factor	1.00
SCAN RANGE	
FOV	25.6
Slice Thickness	1.0
Location per Slab	180
Overlap Locations	0
ACQ TIMING	
Freq	256
Phase	256
Freq DIR	A/P
NEX	1.00
Phase FOV	0.90
Auto Shim	Auto
Phase Correction	No
USER CVS	
User CV23	100.00
MULTI-PHASE	
# of Acquisition	0
Seperate Series	0
Mask Phase	0
Mask Pause	0
ASSET	
Slice Acceleration Factor	1.00
Phase Acceleration Factor	1.00
CONTRAST	
Contrast Yes/No	No