

Protocol: EPILEPSY SETTINGS

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Gradient Mode	ZOOM
Patient Position	Supine	Imaging Mode	2D
Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo
Plane	AXIAL	Imaging Options	MPh, EPI
Series Description	fMRI Resting State 5:00	Acceleration Factor	1.00
SCAN TIMING		SCAN RANGE	
Flip Angle	70	FOV	22.0
TE	Min Full	Slice Thickness	4.0
Number of Echoes	1	Slice Spacing	0.0
TR	2000.0		
Number of Shots	1		
IMAGE ENHANCE		ACQ TIMING	
Filter Choice	None	Freq	64
		Phase	64
		Freq DIR	R/L
		NEX	1.00
		Phase FOV	1.00
		Auto Shim	Auto
		Phase Correction	Yes
GATING/TRIGGER		USER CVS	
Auto Trigger Type	Off	User CV0	1.00
FMRI		MULTI-PHASE	
Initial State	Stimulus	Slice per Location	150
PSD Trigger	Internal	Phase Acquisition Order	Interleaved
Slice Order	Interleaved	Delay after Acquisition	Minimum
View Order	Bottom/Up	# of Acquisition	0
# of Repetitions REST	15	Seperate Series	0
# of Repetitions ACTIVE	15	Delay after Acquisition without AV	0
		Mask Phase	0
		Mask Pause	0
SAT		ASSET	
Tag Type	None	Slice Acceleration Factor	1.00
		Phase Acceleration Factor	1.00
TRICKS		TRACKER	
Pause On/Off	On	Tracker Length	200.0
Auto Subtract	0	Tracker Thickness	20.0
Auto SCIC	Off		
CONTRAST		MULTI-STATION	
Contrast Yes/No	No	Number of Stations	0
		Station Number	0

fMRI Resting State 5:00

fMRI Resting State 5:00