

INFO PAGE		GEOMETRY		CONTRAST	
Total scan duration	04:25.6	Nucleus	H1	Scan type	Imaging
Rel. signal level (%)	100	Coil selection 1	SENSE-Head-32P	Scan mode	3D
Act. TR/TE (ms)	8.2 / 3.7	element selection	selection 1	technique	FFE
ACQ matrix M x P	212 x 172	Coil selection 2	SENSE-Head-32AH	Contrast enhancement	T1
ACQ voxel MPS (mm)	1.00 / 1.00 / 1.00	element selection	selection 1	Acquisition mode	cartesian
REC voxel MPS (mm)	0.95 / 0.96 / 1.00	Dual coil	yes	3D cylindrical mode	no
Scan percentage (%)	100	CLEAR	yes	Fast Imaging mode	TFE
TFE shots	89	body tuned	yes	3D non-selective	no
TFE dur. shot / acq (ms)	1622.9 / 1428.6	FOV AP (mm)	212	+Excitation pulse	default
Min. TI delay	753.2041	RL (mm)	172	shot mode	multishot
Act. WFS (pix) / BW (Hz)	2.252 / 192.9	FH (mm)	150	TFE factor	175
Min. WFS (pix) / Max. BW (Hz)	0.556 / 781.1	Voxel size AP (mm)	1	3D free factor	no
SAR / head	< 6 %	RL (mm)	1	startup echoes	default
Whole body / level	0.0 W/kg / normal	FH (mm)	1	shot interval (ms)	user defined 3000
B1 rms	0.56 uT	Recon voxel size (mm)	1	profile order	linear
PNS / level	60 % / normal	RFOV (%)	81.13207	turbo direction	Z
Sound Pressure Level (dB)	19.06217	Fold-over suppression	no	Echoes	1
MOTION		Slice oversampling	default	partial echo	no
Cardiac synchronization	no	RF select. FOS	no	shifted echo	no
Heart rate > 250 bpm	no	Matrix scan	212	TE	shortest
Respiratory compensation	no	reconstruction	224	Flip angle (deg)	8
Navigator respiratory comp	no	Scan percentage (%)	100	TR	shortest
Flow compensation	no	SENSE	yes	Halfscan	no
fMRI echo stabilisation	no	P reduction (RL)	2	Water-fat shift (pixels)	user defined 2.25
Motion smoothing	no	P os factor	1	Shim	auto
NSA	1	S reduction (FH)	1.2	mDIXON	no
+PROMO	no	k-t BLAST	no	Fat suppression	no
DYN/ANG		Overcontiguous slices	no	Water suppression	no
Angio / Contrast enh.	no	Stacks	1	TFE prepulse	invert
Quantitative flow	no	slices	150	+PP use crusher	no
CENTRA	no	slice thickness (mm)	1	slice selection	no
Manual start	no	slice orientation	transverse	delay	user defined
Dynamic study	no	fold-over direction	RL	To Exc (ms)	900
Arterial Spin labeling	no	fat shift direction	P	PSIR	no
vascular crushing	no	Chunks	1	+tfe switch	no
+Mag reset	no	PlanAlign	no	+MAGIC enable	no
POST/PROC		REST slabs	0	MTC	no
Preparation phases	auto	Catheter tracking	no	+FLEX	no
Interactive F0	no	Interactive positioning	no	T2prep	no
SmartPlan survey	no	External control	no	Research prepulse	no
B0 field map	no	Allow table movement	no	Diffusion mode	no
B1 field map	no	OFFC/ANG		SAR mode	high
MIP/MPR	no	Stacks	1	B1 mode	default
Images	M, no, no, no	Stack Offc. AP (P=+mm)	7.1377	SAR Patient data	auto
Autoview image	M	RL (L=+mm)	-2.926421	PNS mode	low
Calculated images	no, no, no, no	FH (H=+mm)	22.37574	Gradient mode	maximum
Reference tissue	Grey matter	Ang. AP (deg)	1.124372	SoftTone mode	no
Preset window contrast	soft	RL (deg)	0		
Reconstruction mode	real time	FH (deg)	0		
Save raw data	no				
Hardcopy protocol	no				
Ringing filtering	rectangular				
Geometry correction	default				
Elliptical k-space shutter	default				