

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Spin Echo
Coil Configuration	8HRBRAIN	Imaging Options	EPI, DIFF, Asset
Plane	AXIAL		
Series Description	DTI		
SCAN TIMING		SCAN RANGE	
TE	Minimum	FOV	24.0
TR	8500.0	Slice Thickness	2.0
Number of Shots	1	Slice Spacing	0.0
IMAGE ENHANCE		ACQ TIMING	
Filter Choice	None	Freq	128
		Phase	128
		Freq DIR	R/L
		Phase FOV	1.00
		Auto Shim	Auto
		Phase Correction	Yes
GATING/TRIGGER		USER CVS	
Auto Trigger Type	Off	User CV5	1.00
FMRI		MULTI-PHASE	
PSD Trigger	Internal	# of Acquisition	0
Slice Order	Interleaved	Seperate Series	0
View Order	Bottom/Up	Mask Phase	0
# of Repetitions REST	0	Mask Pause	0
# of Repetitions ACTIVE	0		
SAT		DIFFUSION	
Tag Type	None	Optimized TE	Yes
Fat/Water Saturation	Fat	Diffusion Directions	Tensor
		Number of Diffusion Directions	61
		Number of T2 Images	1
		Dual Spin Echo	On
		Diffusion Tensor	No Selection
		Processing Output	
		Recon All Images	Off
TRICKS		CONTRAST	
Pause On/Off	On	Contrast Yes/No	No
Auto Subtract	0		
Auto SCIC	Off		