

\USER\Dr. Behrmann\routine\llan\ep2d_bold_PMU_resting

TA: 8:06 PAT: Off Voxel size: 3.0×3.0×3.0 mm Rel. SNR: 1.00 USER: ep2d_bold_PMU

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	28
Dist. factor	50 %
Position	Isocenter
Orientation	T > C-30.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	2000 ms
TE	30 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize, Elliptical filter
Coil elements	HEA;HEP

Contrast

MTC	Off
Flip angle	73 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	240
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Matrix Coil Mode	Auto (CP)
Distortion Corr.	Off
Unfiltered images	Off
Prescan Normalize	On
Raw filter	<On>
Elliptical filter	On
Mode	Inplane
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

None

System

Body	Off
HEP	On
HEA	On
SP4	Off
SP2	Off
SP8	Off
SP6	Off
SP3	Off
SP1	Off
SP7	Off
SP5	Off

Positioning mode

REF

Table position	H Table
position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Auto Coil Select	Default

Shim mode

Standard

Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	T > C-30.0
Rotation	0.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	125 mm

Physio

1st Signal/Mode None

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active

SIEMENS MAGNETOM Verio syngo MR B17

Meas[16]	Active Meas[17]
	Active Meas[18]
	Active Meas[19]
	Active Meas[20]
	Active Motion
correction	Off Spatial filter
	Off

Sequence

Introduction	On Bandwidth	1628
Hz/Px Free echo spacing	Off	
Echo spacing	0.68 ms	

EPI factor	64	
RF pulse type	Normal	
Gradient mode	Fast	

Physio_Log	On	

\USER\Dr. Behrmann\routine\llan\ep2d_bold_PMU_resting

TA: 8:06 PAT: Off Voxel size: 3.0×3.0×3.0 mm Rel. SNR: 1.00 USER: ep2d_bold_PMU

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	21
Dist. factor	50 %
Position	Isocenter
Orientation	T > C-30.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	1500 ms
TE	30 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize, Elliptical filter
Coil elements	HEA;HEP

Contrast

MTC	Off
Flip angle	73 deg
Fat suppr.	Fat sat.

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	240
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

PAT mode	None
Matrix Coil Mode	Auto (CP)

Distortion Corr.	Off
Unfiltered images	Off
Prescan Normalize	On
Raw filter	<On>
Elliptical filter	On
Mode	Inplane
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

None

System

Body	Off
HEP	On
HEA	On
SP4	Off
SP2	Off
SP8	Off
SP6	Off
SP3	Off
SP1	Off
SP7	Off
SP5	Off

Positioning mode	REF
Table position	H Table
position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	T > C-30.0
Rotation	0.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	125 mm

Physio

1st Signal/Mode None

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active

SIEMENS MAGNETOM Verio syngo MR B17

Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	On
Bandwidth	1628 Hz/Px
Free echo spacing	Off
Echo spacing	0.68 ms

EPI factor	64
RF pulse type	Normal
Gradient mode	Fast

Physio_Log	On

