

PHILIPS INTERA 3T

```

Coil selection = "SENSE-Head-8";
  element selection = "SENSE";
  connection = "d";
Dual coil = "no";
CLEAR = "yes";
  body tuned = "no";
FOV      RL (mm) = 230;
         AP (mm) = 230;
         FH (mm) = 128;
Voxel size  RL (mm) = 3.58999991;
            AP (mm) = 3.59375;
Slice thickness (mm) = 4;
Recon voxel size (mm) = 3.59375;
Fold-over suppression = "no";
Reconstruction matrix = 64;
SENSE = "yes";
  P reduction (AP) = 2.00000024;
  P os factor = 1;
k-t BLAST = "no";
Stacks = 1;
  type = "parallel";
  slices = 32;
  slice gap = "user defined";
    gap (mm) = 0;
  slice orientation = "transverse";
  fold-over direction = "AP";
  fat shift direction = "P";
Stack Offc. AP (P=+mm) = -9.82026482;
           RL (L=+mm) = -2.4874742;
           FH (H=+mm) = 20.1270466;
  Ang. AP (deg) = -0.511259973;
      RL (deg) = 10.8807707;
      FH (deg) = -4.44921541;
Minimum number of packages = 1;
Slice scan order = "ascend";
Large table movement = "no";
PlanAlign = "no";
REST slabs = 0;
Interactive positioning = "no";
Patient position = "head first";
  orientation = "supine";
Scan type = "Imaging";
Scan mode = "MS";
  technique = "FFE";
Contrast enhancement = "no";
Acquisition mode = "cartesian";
Fast Imaging mode = "EPI";
  shot mode = "single-shot";
Echoes = 1;
  partial echo = "no";
  shifted echo = "no";
TE = "user defined";
  (ms) = 33;
Flip angle (deg) = 90;
TR = "shortest";
Halfscan = "no";
Water-fat shift = "minimum";
Shim = "auto";
Fat suppression = "SPIR";
  strength = "strong";
  frequency offset = "default";
Water suppression = "no";
MTC = "no";
Research prepulse = "no";
Diffusion mode = "no";
SAR mode = "high";

```

```

MRI_Protocol_REST.txt
B1 mode = "default";
PNS mode = "moderate";
Gradient mode = "maximum";
SoftTone mode = "no";
Cardiac synchronization = "no";
Respiratory compensation = "no";
Flow compensation = "no";
Temporal slice spacing = "equidistant";
fMRI echo stabilisation = "no";
NSA = 1;
Angio / Contrast enh. = "no";
Quantitative flow = "no";
Manual start = "yes";
Dynamic study = "individual";
    dyn scans = 250;
    recon multiplier = 1;
    dyn scan times = "shortest";
    dummy scans = 4;
    immediate subtraction = "no";
    fast next scan = "no";
    synch. ext. device = "yes";
        start at dyn. = 1;
        interval (dyn) = 1;
    dyn stabilization = "yes";
    prospect. motion corr. = "no";
Keyhole = "no";
Arterial spin labeling = "no";
Preparation phases = "auto";
Manual Offset Freq. = "no";
SENSE ref. scan = "no";
SmartPlan survey = "no";
B0 field map/Dixon = "no";
B1 field map = "no";
MIP/MPR = "no";
Images = " M", (3) " no";
Autoview image = " M";
Calculated images = (4) " no";
Reference tissue = "white matter";
EPI 2D phase correction = "no";
Preset window contrast = "soft";
Reconstruction mode = "real time";
    reuse memory = "no";
Save raw data = "no";
Hardcopy protocol = "no";
Ringing filtering = "default";
Geometry correction = "default";
IF_info_seperator = 0;
Total scan duration = "07:06.7";
Rel. signal level (%) = 100;
Act. TR/TE (ms) = "1667 / 33";
Dyn. scan time = "00:01.7";
ACQ matrix M x P = "64 x 64";
ACQ voxel MPS (mm) = "3.59 / 3.59 / 4.00";
REC voxel MPS (mm) = "3.59 / 3.59 / 4.00";
Scan percentage (%) = 100;
Packages = 1;
Min. slice gap (mm) = 0;
EPI factor = 35;
Act. WFS (pix) / BW (Hz) = "6.701 / 64.8";
BW in EPI freq. dir. (Hz) = "3208.1";
Min. WFS (pix) / Max. BW (Hz) = "6.656 / 65.3";
Min. TR/TE (ms) = "1656 / 9.0";
SAR / head = "< 38 % / 1.2 w/kg";
Whole body / level = "< 0.1 w/kg / normal";
B1 rms [uT] = 1.43506896;
PNS / level = "79 % / normal";
Sound Pressure Level (dB) = 25.5403595;

```