

Resting State Acquisition - UM				
Gradient Echo Reverse Spiral Acquisition (Glover & Law, 2001)				
Voxel Size: 3.438 X 3.438 X 3.0				
Patient Position				
Patient Position:	Supine			
Patient Entry:	Head First			
Coil:	HEAD			
Imaging Parameters				
Plane:	Axial			
Mode:	2D			
Pulse Seq:	Gradient Echo			
Imaging Options:	None			
Scan Timing				
# of Echoes:	1			
TE:	30 ms			
TR:	2000 ms			
Flip Angle:	90			
Acquisition Timing				
Frequency:	64			
Frequency DIR:	A/P			
Auto Center Frequency:	Water			
NEX	1.00			
Autoshim:	On			
Scanning Range				
FOV:	22			
Slice Thickness:	3.0			
Spacing:	0.0			
# of Slices:	40			
User Control Variables				
			Minimum	Maximum
CV 0	# of interleaves (1,2,4,8,16)	1.00	1	128
CV 1	number of temporal frames	300.00	1	16384
CV 2	external trigger (1) or not (0)	0.00	0	1
CV 3	xfer data (script_number) or not (0)	0.00	0	999
CV 4	recon size (64,128,256)	64.00	64	512
CV 5	number of disdaqs	4.00	0	1000