

| PATIENT POSITION | | IMAGING PARAMETERS | |
|-------------------------|-------------|---------------------------|-----------------------|
| Patient Entry | Head First | Imaging Mode | 3D |
| Patient Position | Supine | Pulse Sequence | SPGR |
| Coil Configuration | 8HRBRAIN | Imaging Options | EDR, Fast, IrP, Asset |
| Plane | SAGITTAL | Acceleration Factor | 1.00 |
| Series Description | FSPGR_Anat | | |
| SCAN TIMING | | SCAN RANGE | |
| Flip Angle | 8 | FOV | 25.6 |
| TE | Min Full | Slice Thickness | 1.0 |
| Number of Echoes | 1 | Location per Slab | 176 |
| TI | 600 | Overlap Locations | 0 |
| Receiver Bandwidth | 31.25 | | |
| IMAGE ENHANCE | | ACQ TIMING | |
| Filter Choice | None | Freq | 256 |
| | | Phase | 192 |
| | | Freq DIR | S/I |
| | | NEX | 1.00 |
| | | Phase FOV | 1.00 |
| | | Auto Shim | Auto |
| | | Phase Correction | No |
| GATING/TRIGGER | | USER CVS | |
| Auto Trigger Type | Off | User CV23 | 100.00 |
| FMRI | | MULTI-PHASE | |
| Initial State | Control | # of Acquisition | 0 |
| PSD Trigger | Internal | Seperate Series | 0 |
| Slice Order | Interleaved | Mask Phase | 0 |
| View Order | Bottom/Up | Mask Pause | 0 |
| # of Repetitions REST | 0 | | |
| # of Repetitions ACTIVE | 0 | | |
| SAT | | ASSET | |
| Tag Type | None | Slice Acceleration Factor | 1.00 |
| | | Phase Acceleration Factor | 2.00 |
| TRICKS | | CONTRAST | |
| Pause On/Off | On | Contrast Yes/No | No |
| Auto Subtract | 0 | | |
| Auto SCIC | Off | | |

FSPGR_Anat

FSPGR_Anat