

	Anatomical	Rest
Manufacturer	GE	GE
Model	Discovery	Discovery MR750
Headcoil	8 Chan	8 Chan
Field Strength	3T	3T
Sequence	3D SPGR	EPI
Flip Angle [Deg]	8.0	90
Inversion Time (TI) [ms]	450	-
Echo Time (TE) [ms]	3.136	30
Repetition Time (TR) [ms]	8068.0	2000
Bandwidth per Voxel (Readout) [Hz]	31.25	250
Parallel Acquisition	ASSET x 2	-
Partial Fourier	Off	-
Number of Slices	250	37
Slice Orientation	Saggital	Axial
Slice Phase Encoding Direction	Anterior to Posterior	Anterior to Posterior
Slice Acquisition Order	Interleaved Ascending	Interleaved Ascending
Slice Thickness [mm]	1.0	3.5
Slice Gap [mm]	0	0
Field of View [mm]	250	224
Acquisition Matrix	250x250	64x64
Slice In-Place Resolution [mm <sup>2</sup> ]	1.0x1.0	3.5x3.5
Number of Measurements	-	180
Acquisition Time [min:sec]	5:01	6:04
Fat Suppression	None	Yes
Number of Directions	-	-
Number of B Zeros	-	-
B Value(s) [s/mm <sup>2</sup> ]	-	-
Averages	-	-