

	Anatomical	Rest
Manufacturer	Siemens	Siemens
Model	TrioTim	TrioTim
Headcoil	8 Chan	8 Chan
Field Strength	3T	3T
Sequence	3D MPRAGE	EPI
Flip Angle [Deg]	8.0	80
Inversion Time (TI) [ms]	900	-
Echo Time (TE) [ms]	3.02	30
Repetition Time (TR) [ms]	2600.0	2500
Bandwidth per Voxel (Readout) [Hz]	180	2240
Parallel Acquisition	-	-
Partial Fourier	6/8	-
Number of Slices	176	38
Slice Orientation	-	Axial
Slice Phase Encoding Direction	Anterior to Posterior	Anterior to Posterior
Slice Acquisition Order	Sequential Ascending	Interleaved Ascending
Slice Thickness [mm]	1.0	3.0
Slice Gap [mm]	0.5	0.33
Field of View [mm]	256	216
Acquisition Matrix	256x256	72x72
Slice In-Place Resolution [mm <sup>2</sup> ]	1.0x1.0	3.0x3.0
Number of Measurements	-	184
Acquisition Time [min:sec]	8:19	7:45
Fat Suppression	None	Yes
Number of Directions	-	-
Number of B Zeros	-	-
B Value(s) [s/mm <sup>2</sup> ]	-	-
Averages	-	-